





















# ACETOCHLOR DETAILED FLOWCHART

<p>1.</p>  <p>Remove upper rack from magnetic base. Label test tubes for Standards, Control, and Samples.</p> <table border="1"> <thead> <tr> <th>Tube #</th> <th>Content</th> </tr> </thead> <tbody> <tr> <td>1, 2</td> <td>Diluent/Zero Standard 0 ppb</td> </tr> <tr> <td>3, 4</td> <td>Standard 1, 0.1 ppb</td> </tr> <tr> <td>5, 6</td> <td>Standard 2, 0.5 ppb</td> </tr> <tr> <td>7, 8</td> <td>Standard 3, 2.5 ppb</td> </tr> <tr> <td>9</td> <td>Control</td> </tr> <tr> <td>10</td> <td>Sample 1</td> </tr> <tr> <td>11</td> <td>Sample 2</td> </tr> <tr> <td>12</td> <td>Sample 3</td> </tr> </tbody> </table> <p>Add 200 <math>\mu</math>L of either Standards, Control or Samples to the bottom of each test tube by inserting the pipette tip all the way into the bottom of the tube</p>	Tube #	Content	1, 2	Diluent/Zero Standard 0 ppb	3, 4	Standard 1, 0.1 ppb	5, 6	Standard 2, 0.5 ppb	7, 8	Standard 3, 2.5 ppb	9	Control	10	Sample 1	11	Sample 2	12	Sample 3	<p>6.</p>  <p><b>Do not</b> separate upper rack from lower base. Using a smooth motion, <i>invert</i> the combined rack assembly over a sink and pour out the tube contents; keep inverted and <b>gently blot</b> the test tube rims on several layers of paper toweling.</p>
Tube #	Content																		
1, 2	Diluent/Zero Standard 0 ppb																		
3, 4	Standard 1, 0.1 ppb																		
5, 6	Standard 2, 0.5 ppb																		
7, 8	Standard 3, 2.5 ppb																		
9	Control																		
10	Sample 1																		
11	Sample 2																		
12	Sample 3																		
<p>2.</p>  <p>Add 250 <math>\mu</math>L of Acetochlor Enzyme Conjugate down the inside wall of each tube by aiming the pipet tip 1/4" to 1/2" below the tube rim without touching the rim or tube wall with the pipet tip; deliver liquid gently.</p>	<p>7.</p>  <p>Add 1 mL of Washing Solution down the inside wall of each tube by using the technique described in Box 2. <i>Wait 2 minutes</i>. Using a smooth motion, invert the combined rack assembly over a sink and pour out the tube contents: keep inverted and <b>gently blot</b> the test tube rims on several layers of paper toweling. Repeat this step.</p>																		
<p>3.</p>  <p>Add 500 <math>\mu</math>L of thoroughly mixed Acetochlor Antibody Coupled Magnetic Particles down the inside wall of each tube by using the technique described in Box 2. <i>Vortex</i> for 1 to 2 seconds (at low speed to minimize foaming).</p>	<p>8.</p>  <p>Lift the upper rack (with its tubes) off the magnetic base; add 500 <math>\mu</math>L of Color Reagent down the inside wall of each tube by using the technique described in Box 2. <i>Vortex</i> for 1 to 2 seconds (at low speed to minimize foaming).</p>																		
<p>4.</p>  <p>React 20 minutes at room temperature ( 15° - 30°C).</p>	<p>9.</p>  <p>React for 20 minutes at room temperature (15° - 30° C). During this period, add 1 mL of Washing Solution into a clean tube for use as an instrument blank in Step 10.</p>																		
<p>5.</p>  <p>Combine the upper rack with the magnetic base; press all tubes into base; allow 2 minutes for the particles to separate.</p>	<p>10.</p>  <p>Add 500 <math>\mu</math>L of Stopping Solution down the inside wall of each tube by using the technique previously described. <i>Read</i> results at 450 nm within 15 minutes after adding the Stopping Solution. <i>Multiply</i> results of samples by the appropriate dilution factor (if any).</p> <p><b>[Safety Caution:</b> Stopping Solution contains diluted sulfuric acid.]</p>																		

For Ordering or Technical Assistance Contact:  
 ABRAXIS, LLC 54 Steamwhistle Drive, Warminster, PA 18974  
 Phone: 215-357-3911 Fax: 215-357-5232  
 Web: [www.abraxiskits.com](http://www.abraxiskits.com)

Acetochlor Magnetic Particle Kit Part # 500021, 100 Test

# ACETOCHLOR CONCISE FLOWCHART

<p>1.</p>  <p>Separate the rack.</p> <p>Add 200 <math>\mu</math>L of either Standards, Control or Samples to the bottom of each test tube.</p>	<p>6.</p>  <p>Invert the combined rack.</p> <p>Blot <b>gently</b>.</p>
<p>2.</p>  <p>Add 250 <math>\mu</math>L of Acetochlor Enzyme Conjugate to each test tube.</p>	<p>7.</p>  <p>Add 1 mL of Washing Solution.</p> <p>Wait 2 minutes.</p> <p>Invert the combined rack.</p> <p>Blot <b>gently</b>.</p> <p>Repeat this step.</p>
<p>3.</p>  <p>Add 500 <math>\mu</math>L of mixed Magnetic Particles to each test tube.</p> <p>Vortex.</p>	<p>8.</p>  <p>Separate the rack.</p> <p>Add 500 <math>\mu</math>L of Color Reagent to each test tube.</p> <p>Vortex.</p>
<p>4.</p>  <p>Incubate for 20 minutes.</p>	<p>9.</p>  <p>Incubate for 20 minutes.</p> <p>Prepare blank.</p>
<p>5.</p>  <p>Combine the rack and magnetic base.</p> <p>Seat all tubes.</p> <p>Wait 2 minutes.</p>	<p>10.</p>  <p>Add 500 <math>\mu</math>L of Stopping Solution to each test tube.</p> <p>Read OD 450</p>

For Ordering or Technical Assistance Contact:  
 ABRAXIS, LLC 54 Steamwhistle Drive, Warminster, PA 18974  
 Phone: 215-357-3911 Fax: 215-357-5232  
 Web: [www.abraxiskits.com](http://www.abraxiskits.com)

Acetochlor Magnetic Particle Kit Part # 500021, 100 Test