



Rabbit anti-*E. coli* O104:H4
(Polyclonal Antibody)

Phone: (215) 357-3911
FAX: (215) 357-5232
Email: info@abraxiskits.com
WEB: Abraxiskits.com

Target: *E. coli* O104:H4

Catalog #: 412030

Immunogen: *E. coli* O104:H4

Specificity: reacts with *E. coli* O104:H4

Clonality: Polyclonal

Species: Raised in rabbit


Purity: Protein-A/G Affinity purified

Form: 0.1 mg (as determined by UV absorbance at 280nm) lyophilized, reconstitute in 1x PBS.

Applications: ELISA IMS Beads
Dipsticks LAT
Immuno PCR Labeling

Storage: Store at 4 °C after reconstitution for up to 1 week. Aliquot and freeze at -20 °C to -80 °C for longer term storage. Avoid repeated freeze/thawing.

For Laboratory research only, not for drug, diagnostics or other use.



Cross Reactivity profile of *E. coli* O104 Polyclonal Antibody against 6 STEC and O157 with 10^7 cells coated onto microtiter plates using an indirect ELISA.